



## INSTRUCTIONS BLOOD TEST

*This is a reliable test voor syphilis and HIV*


Please take your time to read the instructions. It shouldn't be difficult, but please follow the instructions carefully. This way, errors between the collection and sending to the lab can be prevented.




1. Put the collection tube, the finger stick, and the spare stick on a clean surface. Open the tube by removing the pink cap.



2. Wash your hands with warm water and soap. If your hands are cold, warm them up in hot water first. Massage your middle finger to make the blood flow to the top.



3. Take the finger stick and remove the cap.



4. Place your hand back down on the table. Place the stick firmly on the inside of the fingertip.



5. You should hear a click, which immediately is followed by the finger prick. Don't remove the stick until you heard the click.



6. Massage with your other hand as much blood as possible to the fingertip and let it drop into the tube. Try to fill the tube with blood as much as possible, **up to the mark of '500'**.



7. Close the cap carefully till you hear a click. Put the tube into the plastic safety bag.

**Close and seal the safety bag only when all tests are inside.**